



MIQE Guidelines

Minimum Information for
Publication of Quantitative Real-
Time PCR Experiments

產品專員 / 陳立德


正茂生物科技有限公司
Genmali Biotechnology Co., Ltd.

TEL : (02) 2731-3369 FAX : (02) 2731-9590
免費諮詢服務及訂貨專線 : 0800-213-029



Central Dogma

Central dogma describes
information flow from
DNA→RNA→protein




www.bio-rad.com/pcr





DNA Structure



Deoxyribonucleic Acid (DNA)

Phosphoric acid backbone

Base pairs

Nitrogenous Bases

Guanine (G) Cytosine (C)

Adenine (A) Thymine (T)

Hydrogen bonds

Nucleotide

Base pair

www.bio-rad.com/pcr






Polymerase Chain Reaction

The polymerase synthesis a new DNA complementary to template in 5' to 3' direction

Primers are complementary to opposite strands of target region but not complementary to any other sequences

EXTENDED PRIMERS



95°C

50-65°C

72°C

72°C

50-65°C



www.bio-rad.com/pcr






PCR – 1st Generation

- Amplify target DNA with end point analysis to distinguish products
 - No relationship between end point and starting target copies

End point analysis

www.bio-rad.com/pcr





What is Real-Time PCR?

Simply - **fluorescent molecules** are used to monitor the reaction while amplification is taking place.

You are able to view this occurring in **real-time** on your instrument.

www.bio-rad.com/pcr






Reality vs. Theory

Amplification is exponential, but the exponential increase is limited:

- A linear increase follows exponential
- Eventually plateaus

Real-Time PCR allows us to 'see' the exponential phase so we can calculate how much we started with.



Theoretical

Real Life


www.bio-rad.com/pcr





Quantitative PCR (qPCR)

Real-time qPCR enables assessment of the reaction after each cycle










Product Amount

Cycle number

Conventional PCR

www.bio-rad.com/pcr



What are the most common detection strategies used for Real-Time PCR?

www.bio-rad.com/pcr

Real-Time PCR

- These fluorescent molecules can be used
 - Non-specific DNA binding dyes
 - SYBR® Green I
 - Ethidium Bromide
 - Specific Hybridization Probes/Primers
 - TaqMan™
 - molecular beacons
 - dual-oligo FRET pairs
 - Scorpions™/Amplifluor™/LUX™

www.bio-rad.com/pcr

DNA binding dyes

Add Master Mix & Sample

Reaction Tube

Denaturation

Annealing

www.bio-rad.com/pcr

DNA binding dyes

Extension

Extension Continued

Apply Excitation Wavelength

Repeat





www.bio-rad.com/pcr

Melt Curve Analysis

Fluorescence decrease as the temperature increase:

- DNA strands start to separate
- SYBR green loses its binding to the DNA
- Fluorescence rapidly decreases

www.bio-rad.com/pcr



AMPLIFICATION

Probe-based assay

- More specific than dye base method
- Multiplex able using different reporter dye
- SNP genotyping application
- More high cost than dye base method

www.bio-rad.com/pcr

AMPLIFICATION

Historically... before qPCR

10 fold! I don't see this as 10 fold. It looks more like 3 fold to me!

- Laborious techniques were used to examine specific targets.
- These methods were Semi-quantitative at best.

www.bio-rad.com/pcr
Clinical Chemistry 55:4 (2009) 613

AMPLIFICATION

After the arrival of qPCR

- Assay development was much easier.
- Sample required was much smaller.
- Results were achieved much faster.
- Data sets grew larger.
- Instruments generated numerical outputs.

www.bio-rad.com/pcr

AMPLIFICATION

The Problem

- Perceived simplicity of running qPCR experiments has resulted in many publications containing seriously misleading results.

www.bio-rad.com/pcr

AMPLIFICATION

QPCR papers

Pubmed

Summary • 20 per page • Sort by Most Recent

Results: 1 to 20 of 113592

1. Promoter Methylation Status Modulates the Expression of Tumor Suppressor (Rb2p130) Gene in Breast Cancer.
Ullah F, Khan T, Ali N, Malik FA, Kayani MA, Shah ST, Saeed M.
AIDS Res 2013 Aug 13;13(8):851-857. doi: 10.1177/1042445113504851. eCollection 2013.
PMID: 23711334
2. An integrated system for identifying the hidden assassins in traditional medicines containing antitubercular agents.
Ali L, Sun W, Wang B, Zhao H, Li Y, Cai S, Xiang L, Zhu Y, Yao H, Song J, Cheng YC, Chen S.
Sci Rep. 2015 Apr 13;5:10111. doi: 10.1038/srep10111.
PMID: 25719549
3. 16S rDNA Gene Sequence-Based Identification of Bacteria in Automatically Incubated Blood Cultures.
Frickenmann H, Denker D, Schwarz NG, Pattn A, Boonen K, Sarpang N, Adi-Sarkodie Y, Hagedewachs E, Wiersma F, von Kries R, Poppert S, Lodenstein U, Mly J, Hagen RM.
AIDS Res 2015 Aug 13;13(8):851-857. doi: 10.1177/1042445113504851. eCollection 2013.
PMID: 25719531
4. Development and Characterization of Probe-Based Real Time Quantitative RT-PCR Assays for Detection and Serotyping of Foot-and-mouth Disease Viruses Circulating in West Eurasia.
Jamaal SM, Rashedi M.
PLoS One 2015 Aug 11;10(8):e0135655. doi: 10.1371/journal.pone.0135655. eCollection 2015.
PMID: 26270102

www.bio-rad.com/pcr

AMPLIFICATION

Questionable Data

AIDS RESEARCH AND HUMAN IMMUNODEFICIENCY VIRUS

Slope : -1.365
Eff = 440.26%

Human Endogenous Retroviral Infection

Slope : -2.276
Eff = 175.02%

FIG. 1. Syncytin-1 DNA copy numbers are increased in brain of MS patients. A standard curve for DNA copy number (correlation coefficient: 0.97; slope: -1.365; intercept: 29.435; $y = -1.365x + 29.435$) was derived (A). The number of viral DNA copies indicated a significant increase in brain of MS patients (9.8 log₁₀) relative to non-MS controls (7.9 log₁₀) (B). A standard curve for GAPDH (correlation coefficient: 0.96; slope: -2.276; intercept: 29.371; $y = -2.276x + 29.371$) was derived (C).

www.bio-rad.com/pcr

AMPLIFICATION

Retracted Paper

RESEARCH ARTICLE

9 SEPTEMBER 2005 VOL 309 SCIENCE www.sciencemag.org

Retraction

WE WISH TO RETRACT OUR RESEARCH ARTICLE "THE MINA OF THE ARABIDOPSIS GENE *FT* MOVES FROM LEAF TO SHOOT APICES AND INDUCES FLOWERING" (1). After the first author (TH) left the Umeå Plant Science Centre for another position, analysis of his original data revealed several anomalies. It is apparent from these files that data from the real-time RT-PCR were analysed incorrectly. Certain data points were removed, while other data points were given increased weight in the statistical analysis. When all the primary real-time RT-PCR data are subjected to correct statistical analysis, most of the reported significant differences between time points disappear. Because of this, we are retracting the paper in its entirety.

In new experiments, we have reproduced the floral induction caused by a heat-shock induction of *FT* in a single leaf, but we have failed to detect movement of the transgenic *FT* mRNA from leaf to shoot apices. We therefore retract the conclusion that *FT* mRNA is part of the floral inductive signal moving from leaf to shoot apices.

We deeply regret any scientific misconceptions that have resulted from the publication of these data.

The first author of the paper (T.H.) strongly objects to the retraction of the paper and has therefore declined to be an author of the retraction.

Our related Science Report on the *CUT1* regulatory module in trees (2) is not affected by this retraction. In this paper, T.H. was involved in the construction and analysis of the *FT::GUS* experiments reported in the Supporting Online Material. These data have been reanalysed and found to be correctly reported.

HELVIG BOHLERUS¹, SVEN ERICSSON², FRANÇOIS PARCY³, OVE NILSSON⁴

¹Umeå Plant Science Centre, Department of Forest Genetics and Plant Physiology, Swedish University of Agricultural Sciences, SE-90183 Jönköping, Sweden; ²Laboratoire de Physiologie Cellulaire et Moléculaire, Département de Biologie, Université de Sherbrooke, 3301 Avenue de la Biologie, Sherbrooke, Québec J1K 2R1, Canada; ³Centre National de la Recherche Scientifique (CNRS), Commissariat à l'Energie Atomique (CEA), Institut National de la Recherche Agronomique (INRA), Université Joseph Fourier, 17 rue des Bains, Saint-Etienne 42023, France; ⁴INRA, UR1213, 33077 Villenave d'Ornon, France

References
1. T. Huang, H. Böhlerius, S. Ericsson, J. Parcy, O. Nilsson, Science 309, 549-550 (2005).
2. H. Böhlerius et al., Science 314, 1461-1464 (2006).

www.bio-rad.com/pcr

AMPLIFICATION

Autism and MMR Vaccine

Amplification - f gene 11-04-01

Abn

Cycle

Complex

Viewer: ΔRn (B)

Reporter: FAM

Threshold Cycle Calculation

Use Threshold: 0.01 (Suggest)

Multi-Δ Sliders: 10.0 (Suggest)

Omit Threshold: 2.0



Baseline

Start: 5 Stop: 15

Update Calculations

	CT	Std Dev
FAM - F4	34.744	0.001
FAM - H2	40.000	0.001

www.bio-rad.com/pcr



What can be done?

www.bio-rad.com/pcr

BIO-RAD

- ### MIBBI, RDML and MIQE
- Promoting coherent minimum reporting guidelines for biological and biomedical investigations: the MIBBI project. *Nature biotechnology*, volume 26, number 8, August 2008, 889-896.
 - RDML: structured language and reporting guidelines for real-time quantitative PCR data. Lefever S, Hellemans J, Pattyn F, Przybylski DR, Taylor C, Geurts R, Untergasser A, Vandesompele J; on behalf of the RDML consortium. *Nucleic Acids Res*. 2009 Feb 17.
 - The MIQE Guidelines: Minimum Information for Publication of Quantitative Real-Time PCR Experiments. Bustin SA, Benes V, Garson JA, Hellemans J, Huggett J, Kubista M, Mueller R, Nolan T, Pfaffl MW, Shipley GL, Vandesompele J, Wittwer CT. *Clin Chem*. 2009 Feb 26.
- www.bio-rad.com/pcr
- BIO-RAD

A new beginning

Clinical Chemistry 55:4
611–622 (2009)

Special Report

The MIQE Guidelines: Minimum Information for Publication of Quantitative Real-Time PCR Experiments

Stephen A. Bustin^{1*}, Vladimir Benes², Jeremy A. Garson^{3,4}, Jan Hellemans⁵, Jim Huggett⁶,
Mikael Kubista^{7,8}, Reinhold Mueller⁹, Tania Nolan¹⁰, Michael W. Pfaffl¹¹, Gregory L. Shipley¹²,
Jo Vandesompele⁵ and Carl T. Wittwar^{13,14}

**57 Essential information
(E) must be submitted with the manuscript**

**28 Desirable information
(D) should be submitted if available**

www.bio-rad.com/pcr

What is MIQE? It's a Checklist

- qPCR community driven guidelines for essential and desired information in literature;
 - Experimental Design
 - Sample Information
 - Nucleic Acid Extraction
 - Reverse Transcription
 - qPCR Target Information
 - qPCR Oligonucleotides
 - qPCR Protocol
 - qPCR Validation
 - Data Analysis

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr

Experimental Design

Item to check	Importance
Experimental design	
Definition of experimental and control groups	E
Number within each group	E
Assay carried out by the core or investigator's laboratory?	D
Acknowledgment of authors' contributions	D

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr

Experimental Replicates

Biological Replicates

Experiment: A, B, C
Control: A, B, C

RT-qPCR Samples

Technical replicates

Gene of Interest: 1 A, 2 A, 3 A, 1 B, 2 B, 3 B, 1 C, 2 C, 3 C
Reference Gene: 1 A, 2 A, 3 A, 1 B, 2 B, 3 B, 1 C, 2 C, 3 C

Gene of Interest: A, A, A, B, B, B, C, C, C
Reference Gene: A, A, A, B, B, B, C, C, C

Gene of Interest: [] [] []
Reference Gene: [] [] []

NTC

www.bio-rad.com/pcr

AMPLIFICATION




Sample Information

	Item to check	Importance
Sample	Description	E
	Volume/mass of sample processed	D
	Microdissection or macrodissection	E
	Processing procedure	E
	If frozen, how and how quickly?	E
	If fixed, with what and how quickly?	E
	Sample storage conditions and duration (especially for FFPE ^b samples)	E

Clinical Chemistry 55:4 (2009) 613


www.bio-rad.com/pcr


AMPLIFICATION




Sample Information – Acquisition

- Source (soil, blood, urine, plant tissues ...)
- Sampling techniques (biopsy materials, single-cell sampling, laser microdissection...)
- Description (percentage of the biopsy made up of tumor cells)
- Mass/Volume
- Homogeneity







Clinical Chemistry 55:4 (2009) 613


www.bio-rad.com/pcr


AMPLIFICATION




Sample Information – Processing and Storage

- Processing procedure
- Storage (-80C vs -20C; duration...)
- Prep time to extraction (-80C)
- Sample degradation







Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr


AMPLIFICATION




Nucleic Acid Extraction

	Item to check	Importance
Nucleic acid extraction	Procedure and/or instrumentation	E
	Name of kit and details of any modifications	E
	Source of additional reagents used	D
	Details of DNase or RNase treatment	E
	Contamination assessment (DNA or RNA)	E
	Nucleic acid quantification	E
	Instrument and method	E
	Purity (A_{260}/A_{280})	D
	Yield	D
	RNA integrity: method/instrument	E
	RIN/RQ1 or C_t of 3' and 5' transcripts	E
	Electrophoresis traces	D
	Inhibition testing (C_q dilutions, spike, or other)	E

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr


AMPLIFICATION

Sample Extraction

- DNA / RNA
 - Aquapure™ Genomic DNA Isolation Kit
 - Aurum™ Total RNA Kits
- Contaminants
 - Starches
 - Lipids
 - Metals
- Extraction contaminants
 - Phenol chloroform
 - Salts
- Sample degradation

www.bio-rad.com/pcr

AMPLIFICATION

www.bio-rad.com/pcr

AMPLIFICATION

Inhibitors copurified

BMC Research Notes

Short Report

Differential susceptibility of PCR reactions to inhibitors: an important and unrecognised phenomenon

Jim F Huggett¹, Tanya Novak¹, Jeremy A Garson², Clare Green¹, Stephen D Morris-Jones¹, Robert F Miller⁴ and Alimuddin Zumla¹

Address: ¹Centre for Infectious Diseases and International Health, Wondolop Institute for Medical Sciences, 41 Cleveland Street, University College London, London, W1T 4JF, UK; ²Centre for Virology, Wondolop Institute for Medical Sciences, University College London, London, W1T 4JF, UK; ³Department of Microbiology, University College London, Wondolop Institute for Medical Sciences, London, W1T 4JF, UK and ⁴Research Department of Infection and Population Health, Division of Population Health, University College London, London, W1T 4JF, UK

Email: Jim F Huggett^{*} - j.f.huggett@ucl.ac.uk; Tanya Novak - t.novak@ucl.ac.uk; Jeremy A Garson - j.garson@ucl.ac.uk; Clare Green - c.green@ucl.ac.uk; Stephen D Morris-Jones - s.d.morris-jones@ucl.ac.uk; Robert F Miller - r.f.miller@ucl.ac.uk; Alimuddin Zumla - a.zumla@ucl.ac.uk

^{*} Corresponding author

Published: 28 August 2008
BMC Research Notes 2008, 1:72 doi:10.1186/1756-1556-1-72

Received: 2 August 2008
Accepted: 28 August 2008

Results: When copurified inhibitors were assessed using two different PCR reactions one reaction appeared to be inhibited whilst the other was not. Further experiments using various concentrations of unextracted urine to inhibit six different PCR reactions revealed that susceptibility to inhibition was highly variable between reactions. Similar results were obtained using EDTA as the PCR inhibitor. We could find no obvious explanation why one reaction should be more susceptible to inhibition than another, although a possible association with amplicon GC content was noted.

www.bio-rad.com/pcr

AMPLIFICATION

Effects of EDTA

Available online at www.sciencedirect.com

SCIENCE @ DIRECT[®]

Analytical Biochemistry 351 (2006) 209–210

Notes & Tips

PUD: A quantitative PCR assay for the detection of inhibitors in nucleic acid preparations

Tania Nolan^a, Rebecca E. Hands^b, William Ogunkolade^b, Stephen A. Bustin^{b,*}

^a School of Cell and Molecular Science, Barts and The London, Queen Mary's School of Medicine and Dentistry, University of London, London E1 1BB, UK

^b Department of Cell and Molecular Science, Barts and The London, Queen Mary's School of Medicine and Dentistry, University of London, London E1 1BB, UK

Received 3 November 2005
Available online 20 February 2006

www.bio-rad.com/pcr

Template DNA Preparation

- Genomic DNA
 - Cut with restriction enzyme that does not cut within amplicon
 - Boil DNA for 10 min and then onto ice
- Plasmid DNA
 - If it doesn't work, linearize plasmid with restriction enzyme that does not cut within amplicon
- cDNA
 - Prior to RT: Treat RNA with RNase-free DNase
 - After RT: use enzyme that has RNaseH activity to digest away RNA from RNA:DNA hybrid

www.bio-rad.com/pcr

Sample and Template Preparation

- Extract and analyze RNA
- Careful quantification is necessary
 - RiboGreen Assay - Quantification
 - NanoDrop – Quantification and purity
 - Experion™ - Quality and quantification

www.bio-rad.com/pcr

Analysis of RNA purity and integrity

Samples 6,7,8 are highly degraded.

Sample	Conc (ng/ul)	A260/280*	A260/230*
Control-no heat	115	1.90	2.44
3 min @ 90C	114	1.93	2.40
5 min @ 90C	115	2.06	2.37
10 min @ 90C	115	2.03	2.37
15 min @ 90C	116	2.02	2.31
1 hr @ 90C	109	1.99	2.18
2hr @ 90C	117	2.00	2.32
4 hr @ 90C	118	1.89	2.23

*Generally accepted ratios (A260/280 and A260/230) for good quality RNA are > 1.8.

www.bio-rad.com/pcr

Total RNA analysis



28S →
18S →

bpp
5000
4000
3000
2000
1500
1000
700
500
300

Lane 1: 100 bp DNA ladder
Lane 2: acceptable quality
Lane 3: heavily degraded total RNA
Lane 4: slightly degraded and DNA- contaminated sample
Lane 5: 1 Kb DNA ladder

1% Agarose gel electrophoresis

www.bio-rad.com/pcr



Reverse Transcription

Item to check	Importance
Reverse transcription	
Complete reaction conditions	F
Amount of RNA and reaction volume	E
Priming oligonucleotide (if using GSP) and concentration	E
Reverse transcriptase and concentration	E
Temperature and time	E
Manufacturer of reagents and catalogue numbers	D
C _q s with and without reverse transcription	D ^c
Storage conditions of cDNA	D

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr

BIO-RAD

RT protocols:

RNA was reversed transcribed,

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

Papers in Press. Published October 31, 2014 as doi:10.1373/clinchem.2014.230615
The latest version is at <http://hwmaint.clinchem.org/cgi/doi/10.1373/clinchem.2014.230615>

Clinical Chemistry 61:1
000–000 (2015)

Molecular Diagnostics and Genetics

Variability of the Reverse Transcription Step: Practical Implications


Stephen Bustin,^{1*} Harvinder S. Dhillon,¹ Sara Kirvell,¹ Christina Greenwood,¹ Michael Parker,¹
Gregory L. Shipley,² and Tania Nolan¹

BACKGROUND: The reverse transcription (RT) of RNA to cDNA is a necessary first step for numerous research and molecular diagnostic applications. Although RT efficiency is known to be variable, little attention has been paid to the practical implications of that variability.

RESULTS: RT efficiency is enzyme, sample, RNA concentration, and assay dependent and can lead to variable correlation between mRNAs from the same sample. This translates into relative mRNA expression levels that generally vary between 2- and 3-fold, although higher levels are also observed.

www.bio-rad.com/pcr

BIO-RAD



AMPLIFICATION

Reverse transcription Strategies

- Oligo dT
- Random primers
- Oligo dT and Random primers
- Gene specific primers
- RNase H
- RNase inhibitors

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

RNase is STRONG





Plasmid extracted with traditional 3-steps and eluted with sterile ddH₂O.

NO RNase

Contain RNase

www.bio-rad.com/pcr

BIO-RAD



AMPLIFICATION

qPCR Target Information

Item to check	Importance
qPCR target information	
Gene symbol	E
Sequence accession number	E
Location of amplicon	D
Amplicon length	E
In silico specificity screen (BLAST, and so on)	E
Pseudogenes, retropseudogenes, or other homologs?	D
Sequence alignment	D
Secondary structure analysis of amplicon	D
Location of each primer by exon or intron (if applicable)	E
What splice variants are targeted?	E

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr

AMPLIFICATION

NCBI

- Free resource
- Blast sequence
- Stack sequence

www.bio-rad.com/pcr

AMPLIFICATION

Blast analysis

Distribution of 120 Blast Hits on the Query Sequence

<http://www.ncbi.nlm.nih.gov/BLAST/Blast.cgi>

www.bio-rad.com/pcr

AMPLIFICATION

Amplicon Secondary Structures

<http://mfold.rna.albany.edu/?q=mfold>

- Bad location for primers
- Good location for primers

www.bio-rad.com/pcr

AMPLIFICATION

Primer A: $\eta = 66.3\%$

Primer B: $\eta = 95.8\%$

Primer A

PCR Base Line Subtracted

10 100 1000 10000

0 2 4 6 8 10 12 14 16 18 20 22 24 26 28 30 32 34 36 38 40 42 44 46 48 50

Time

Primer B

10 100 1000 10000

0 2 4 6 8 10 12 14 16 18 20 22 24 26 28 30 32 34 36 38 40 42 44 46 48 50

Cycle

Forward Primer

Reverse Primer A

Reverse Primer B

110

1

www.bio-rad.com/pcr

BIO-RAD



AMPLIFICATION



qPCR Oligonucleotides

Item to check	Importance
qPCR oligonucleotides	
Primer sequences	E
RTPrimerDB identification number	D
Probe sequences	D ^d
Location and identity of any modifications	E
Manufacturer of oligonucleotides	D
Purification method	D

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr



[illegible]

AMPLIFICATION

qPCR Protocol

qPCR protocol	Item to check	Importance
Complete reaction conditions		E
Reaction volume and amount of cDNA/DNA		E
Primer, (probe), Mg ²⁺ , and dNTP concentrations		E
Polymerase identity and concentration		E
Buffer/kit identity and manufacturer		E
Exact chemical composition of the buffer		D
Additives (SYBR Green I, DMSO, and so forth)		E
Manufacturer of plates/tubes and catalog number		D
Complete thermocycling parameters		E
Reaction setup (manual/robotic)		D
Manufacturer of qPCR instrument		E

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

Vol 465 24 June 2010 doi:10.1038/nature09144

nature

ARTICLES

A coding-independent function of gene and pseudogene mRNAs regulates tumour biology

Laura Poliseno¹*, Leonardo Salmena¹*, Jiangwen Zhang², Brett Carver³, William J. Haveman¹ & Pier Paolo Pandolfi¹

The canonical role of messenger RNA (mRNA) is to deliver protein-coding information to sites of protein synthesis. However, given that microRNAs bind to RNAs, we hypothesized that RNAs could possess a regulatory role that relies on their ability to compete for microRNA binding, independently of their protein-coding function. As a model for the protein-coding-independent role of RNAs, we describe the functional relationship between the mRNAs produced by the *PTEN* tumour suppressor gene and its pseudogene *PTENP1* and the critical consequences of this interaction. We find that *PTENP1* is biologically active as it can regulate cellular levels of *PTEN* and exert a growth-suppressive role. We also show that the *PTENP1* locus is selectively lost in human cancer. We extended our analysis to other cancer-related genes that possess pseudogenes, such as oncogenic *KRAS*. We also demonstrate that the transcripts of protein-coding genes such as *PTEN* are biologically active. These findings attribute a novel biological role to expressed pseudogenes, as they can regulate coding gene expression, and reveal a non-coding function for mRNAs.

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

d

e

f

METHODS SUMMARY

Cell lines were cultured under standard conditions. miRNA overexpression was obtained by transient transfection (si-miRNAs). *PTENP1/PTEN/KRASIP* 3' UTR overexpression was achieved by transient transfection (pCMV expression vectors) or stable infection with MSCV-PIG retroviral constructs. miRNA-target interaction was measured by luciferase reporter assay. *PTENP1*, *PTEN*, *KRAS*, *KRASIP* and mRNA expression level was detected by real-time PCR. Proliferation, foci and soft agar assays were performed according to standard protocols.

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

Recipe, Protocol, Consumable, and Instrument

Data Analysis - Gene Expression Multiplex - Time Course 1.pcrd

	1 RXN	# RXNs
TE (1X)	2.8	252.0
2X IQ PowerMix (→1X)	10	900.0
10uM IL185101F (→0.3 uM)	0.6	54.0
10uM IL18246R (→0.3 uM)	0.6	54.0
10uM FAM-IL18 BLine (→0.2uM)	0.4	36.0
10uM GAPDH 169F (→0.3 uM)	0.6	54.0
10uM GAPDH 290R (→0.3 uM)	0.6	54.0
100uM HEX-GAPDHBLine (→0.2uM)	0.4	36.0
	1440.0	Total
Template	4	

www.bio-rad.com/pcr

BIO-RAD

19

AMPLIFICATION

The DNA binding dye reagents of Bio-Rad

Reagents

- iQ SYBR Green Supermix
- Ssofast EvaGreen Supermix
- iTaq™ Universal SYBR® Green One-Step Kit

Key Features

- Sensitive SYBR Green or EvaGreen
- Simplified reaction set up with ready-to-use 2x reaction mix
- iTaq hot-start DNA polymerase, hot-start in 30 sec or 3 min at 95°C
- Linear results over 7 orders of magnitude

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

The probes based reagents of Bio-Rad

Reagents

- iTaq universal supermix for Probe
- iTaq™ Universal probe One-Step Kit
- iQ Multiplex Powermix

Key Features

- Simplified reaction set up with ready-to-use 2x reaction mix
- hot-start DNA polymerase, hot-start in 30 sec or 3 min at 95°C
- Linear results over 7 orders of magnitude
- Detection of up to 2 or 5 targets, where one target can be as much as 10⁶-fold different in input quantity

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

3rd Generation Reagents:

SsoFast™ Series Supermix

- High Speed
- High Dynamic Range
- High Tolerance

www.bio-rad.com/pcr

BIO-RAD

AMPLIFICATION

3rd Generation Reagents

SsoFast EvaGreen Supermix

Sso7d from *Sulfolobus solfataricus*





- 7kD, 63 aa.
- Thermostable (T_m >90° C)
- No sequence preference
- Binds to dsDNA (3-6 bp/protein molecule)
- Monomeric





Fusing the double-stranded DNA binding protein Sso7d to Proof gives it a powerful sliding grip on the replicated DNA.


- Minimal inhibition of PCR by use of EvaGreen
- Higher activity
- Tolerant to PCR inhibitors

www.bio-rad.com/pcr

BIO-RAD







Inhibitor Concentration Not Having Any Significant Effect on qPCR Reaction

Inhibitor	Max Tested	Dilution Series	EVA	iQ SYBR	ABI	INV	Euro
Ethanol	25 %	2 fold	12.5 %	6.25 %	3.125 %	3.125 %	3.125 %
Isopropanol	17.5 %	2 fold	8.75 %	4.375 %	2.18 %	4.375 %	4.375 %
KCl	500 mM	2 fold	62.5 mM	31.25 mM	125 mM	62.5 mM	31.25 mM
NaAc	250 mM	2 fold	62.5 mM	62.5 mM	62.5 mM	62.5 mM	15.63 mM
Phenol/Chloroform	6.25 %	2 fold	1.53 %	0.781 %	0.781 %	0.781 %	0.391 %
Blood Serum	10 %	4 fold	2.5 %	0.0098 %	0.0098 %	0.0098 %	0.0098 %
Total Blood (+hep)	10 %	4 fold	0.156 %	0.0024 %	0.0098 %	0.0024 %	0.0024 %
Heparin	0.0188 U	2 fold	0.0188 U	0.00469 U	0.00059 U	0.00023 U	0.00012 U
Soil	25 %	4 fold	6.25 %	1.56 %	1.56 %	6.25 %	0.39 %
Coffee	25 %	4 fold	6.25 %	1.56 %	1.56 %	1.56 %	1.56 %

www.bio-rad.com/pcr

BIO-RAD

qPCR Validation

Item to check	Importance
qPCR validation	
Evidence of optimization (from gradients)	D
Specificity (gel, sequence, melt, or digest)	E
For SYBR Green I, C_q of the NTC	E
Calibration curves with slope and y intercept	E
PCR efficiency calculated from slope	E
CLs for PCR efficiency or SE	D
r^2 of calibration curve	E
Linear dynamic range	E
C_q variation at LOD	E
CLs throughout range	D
Evidence for LOD	E
If multiplex, efficiency and LOD of each assay	E

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr





BIO-RAD

Optimization

How to choose temperature?

www.bio-rad.com/pcr

BIO-RAD



AMPLIFICATION




Data Analysis

Item to check	Importance
Data analysis	
qPCR analysis program (source, version)	E
Method of C_q determination	E
Outlier identification and disposition	E
Results for NTCs	E
Justification of number and choice of reference genes	E
Description of normalization method	E
Number and concordance of biological replicates	D
Number and stage (reverse transcription or qPCR) of technical replicates	E
Repeatability (intraassay variation)	E
Reproducibility (interassay variation, CV)	D
Power analysis	D
Statistical methods for results significance	E
Software (source, version)	E
C_q or raw data submission with RDML	D

Clinical Chemistry 55:4 (2009) 613

www.bio-rad.com/pcr


AMPLIFICATION




qPCR Program/Software




CFX Manager

www.bio-rad.com/pcr


AMPLIFICATION




Absolute Quantization




www.bio-rad.com/pcr


AMPLIFICATION




Relative Quantization



www.bio-rad.com/pcr


Gene Expression Normalization - Comparative C_q

- Relative Quantity (ΔC_q)
 - Not normalized
 - Normalization accomplished via equal loading of samples
 - Post analysis normalization
- Normalized Expression ($\Delta\Delta C_q$)
 - Accounts for loading differences
 - Usually normalize to reference gene
 - Relative quantity of C_q is normalized by the relative quantity of the reference genes

www.bio-rad.com/pcr

BIO-RAD

Relative Quantification - ΔC_q

	C _q
Sample :	22
Control :	24
Delta Ct:	22 - 24 = - 2

Fold induction = $2^{-(-2)} = 4$

www.bio-rad.com/pcr

BIO-RAD

Comparative C_q Method ($2^{-\Delta\Delta C_q}$)

- ❖ Normalize C_q of target gene to that of reference gene
- ❖ Normalize ΔC_q of sample to ΔC_q of control sample

	Reference	C _q
Tissue #1: (sample)	21	22
Tissue #2: (control)	20	24
1 st Delta	Delta C _q #1:	22-21 = 1
	Delta C _q #2:	24-20 = 4
2 nd Delta	Delta C _q :	1 - 4 = - 3
	Fold induction =	$2^{-(-3)} = 8$

www.bio-rad.com/pcr

BIO-RAD

Relative Quantification –Pfaffl Modification

	Reference	Target C _q
Tissue #1: (sample)	21	22
Tissue #2: (control)	20	24
(From Standard curve) Efficiency:	90% = 1.9	100% = 2
Delta C _q :	20-21 = -1	24-22 = 2

Fold induction = $\frac{2^{\frac{\text{deltaCq}_{\text{target GOI}} (24-22 = 2)}{1.9_{\text{reference}}}}}{2^{\frac{\text{deltaCq}_{\text{reference}} (20-21 = -1)}{0.53}}} = \frac{4}{0.53} = 7.5$

www.bio-rad.com/pcr

BIO-RAD

Gene Expression Normalization - Normalized Expression

- $\Delta\Delta Cq$
 - Assume 100% efficiency
 - Only one Ref Gene
- Pfaffl Modification
 - Accounts for efficiency differences
 - Only one Ref Gene
- Vandesompele Method
 - Accounts for efficiency differences
 - Allows multiple reference genes for normalization

Simple

↓

Complex

www.bio-rad.com/pcr

Allelic Discrimination

www.bio-rad.com/pcr

Protein Folding analysis

Sample : 2.5 µg of purified protein in 25 µl tested with 2 µl of Sypro Orange (25µl / dupli)

Protocol: melting curve ranging from 15°C to 100°C with increment of 0.5°C

Estimated run time : 1h 35min

Raw data

Melt Peak

Calculation of the negative first derivative

Please change the peak type to "Negative" here!

Plate map

Value of the Tm (inflection point of the melting transition)

www.bio-rad.com/pcr

HRM mutation analysis

Lung cancer screen in KRAS mutation on CFX

SsoFast EvaGreen Supermix Protocol

98°C – 2 mins
(40 cycles)

98°C – 5 sec

61°C – 10 sec
(Melt)

95°C – 30 sec

75°C – 5 sec (+0.2°C/cycle to 95°C)

Temperature Shifted Difference Curve

Normalized Temperature Shifted Melt Curve

Well	Content	Sample	Cluster	Percent Confidence
A06	Unknown	A549-92bp	Cluster 1	99.7
B06	Unknown	A549-92bp	Cluster 1	99.7
C06	Unknown	HCT116_99bp	Cluster 2	96.8
D06	Unknown	HCT116_99bp	Cluster 2	97.2
E06	Unknown	48T-92bp	Cluster 3	96.5
F06	Unknown	48T-92bp	Cluster 3	97.3
G06	NTC	NTC	<Excluded>	N/A
H06	NTC	NTC	<Excluded>	N/A

www.bio-rad.com/pcr

The RDML International Consortium

- Key developer group, a member community and supporters
- **41 supporters and members from 20 different countries**

22 academic
Austria, Switzerland, France
Spain, Israel, Canada, ...
Sweden (2)
Germany (3)
UK (3)
USA (5)

19 companies
USA (6)
UK (2)
Belgium (2)
Portugal, Austria
Switzerland, Israel,
Australia, Korea, ...

Goal: develop and maintain the RDML data exchange format

www.bio-rad.com/pcr

Conforming to MIQE: an example

- The Bio-Rad technote 5859 provides an excellent summary of a MIQE compliant experiment in practice, from start to finish

amplification

A Practical Approach to RT-qPCR, Publishing Data That Conforms to the MIQE Guidelines

for publication of quantitative real-time PCR experiments (MIQE) guidelines have been recently published (Bustin et al. 2009). This has been followed by the development of an XML-based real-time PCR data markup language (PCML) for the consistent reporting of real-time PCR data by the RDML consortium (Lafleur et al. 2009). This consortium is active in the development of appropriate and standardized terminology, guidelines on minimum information for biological and biomedical investigations, and a flexible and universal data file structure with tools to create, process, and validate HUML files. All relevant information about the RDML project is available at www.rdml.org/.

www.bio-rad.com/pcr

Complete solutions

Gene Modulation
RNAi
siRNA
CRISPR

RNA Purification
Aurum Total RNA Kits
siLentFect
Transfectin
Gene Pulser Xcell

RNA Analysis (Quality and Yield)
Experion System
SmartSpec Plus

Reverse Transcription and Amplification
iScript
iQ SYBR T100
SsoFast reagents
CFX Series

Profiling & Quantification
qBASE software

www.bio-rad.com/pcr

TECHNOLOGY FEATURE

Setting up for higher throughput	209
Physical limitations	209
Cycling	211
Data standards	212
Box 1: Melting shifts show genetic variation	208

qPCR: quicker and easier but don't be sloppy

Monya Baker

Gene profiling using quantitative PCR is becoming higher throughput, but researchers must be careful in gathering their data.

NATURE METHODS | VOL.8 NO.3 | MARCH 2011 | 207

www.bio-rad.com/pcr

